


ICARST 2017

1st International Conference on Applications
of Radiation Science and Technology

A Gamma Ray CT Computed Tomography for Investigating the Wood Structure

A. Saadaoui, I. R. Alami, I. Marrhich, A. Ezzahri, S. Mimount,
A. Ouardi, F. Oujebbour, R. Outayad

The National Centre for Nuclear Energy, Science and Technology
CNESTEN, PB 1382, RP 10001 Rabat - MOROCCO



OUTLINE

- Introduction
- Tomography principle
- Materials & system used
- Experiments
- Results
- Conclusion



ICARST 2017

1st International Conference on Applications
of Radiation Science and Technology





Objectives

- Preliminary study of a cork oak trees of the Maamoura Forest
- The difference in densities
- Attempt for the Datation by counting number of the rings of growth

MOROCCO

Maamoura Forest is located near to Rabat



CNESTEN is inside the Maamoura Forest

Maâmoura Forest

- The Maâmora forest is among the largest s cork oak forests in the world and is the lung of the urban area of Kénitra, Rabat-Salé.
- Approximatively 134 thousand hectares of the forest is filled with cork trees.
- More than 50% of Moroccan's cork trees
- Pine groves, eucalyptus trees and wild pear trees are scattered across the area.




Techniques used in CNESTEN

- Dendrochronolgy
- Dendrochronology is used for datation aims
- Recently another technique is used by studying stable isotopes of carbon, oxygen, and/or hydrogen of tree rings to identify the various climatological, ecological, or hydrological contexts.
- Destructive technique

Tomography


- The tomography is a non-destructive imaging technique.
- This technique gives us an cross sectional images « « slices » ».
- The origin of this technique is a Radon Mathematical theorem



$$R[f(x,y)] = \int_0^{\pi} p(u,\theta) d\theta$$

Tomography

- When gamma rays are directed onto an object, some of the photons interact with the particles of the matter and their energy can be absorbed or scattered.
- Other photons travel completely through the object without interacting
- The number of photons transmitted through an object is depending on the thickness, density and atomic number of the material, and the energy of the individual photons.
- Number of detected photons depends also on a detector efficiency & collimation





$$I = I_0 e^{(-\mu x)}$$

Tomography principle

Example of raw data (sinogram) showing a list of measurements (Y-axis) and corresponding values (X-axis).

Y-axis (Measurements)	X-axis (Values)
64	
128	
21975	21712 2159
22139	21736 2177
21883	21720 2207
21695	21784 2166
21908	21921 2199
21930	21734 2190
21785	21687 2163
21721	21679 2156



System used & materials

- GORBIT system – Electronic & mechanical solved problems
- French NaI detectors (TTL output)
- AM-241 & CO-60 sources (5mm collimation)
- iGorbit software for Reconstruction
- Matlab for image treatments
- Cross section of an oak tree



Case studied

- Cork oak tree cross section



Height = 6cm

Experiments

CO-60 (1,17 & 1,32 MEV)

- 8 mCi

AMERI-241 (59 KEV)

- 500 mCi


	Rotation Step °	Translation Step (mm)	Sinogram Dim	Time
1 st	11.25 °	2.5 mm		
2 nd	11.25°	1.25 mm		
3 rd	2.8125°	2.5 mm		

Image Reconstruction


- iGorbit sotware
- FBP method
- Cosine Filter

CO


16 projections/ 2.5 mm






16 projections/ 1.25 mm




64 projections/ 2.5 mm



AM



AM vs CO



Am

Max=11240
Min=260


Ratio=43

Co


Max=2845
Min=878

Ratio=3



Results Traitment





Results Traitment



Edge detection
In Matlab



Matlab Imaging Tools



Thickness obtained

Real scale



Relative scale

	Real values	Measured by Matlab
D1	21 cm	161.74
D2	24 cm	179.53
E	7 mm	6.25

23.28 cm



8 mm

Density Comparison



High density

Density Comparison



High densities

Another image is obtained



- X ray source
- 70 KV
- 3mA
- 30s



Conclusion

- Am is better for studying the wood
- We can measure the thickness of cork oak trees by CT
- We can observe the density variation
- Datation with tomography?? More resolution

Thank you



saadaoui@cnesten.org.ma
saadaoui@gmail.com
CNESTEN/DAI

The National Centre for Nuclear Energy, Science and Technology
CNESTEN, PB 1382, RP 10001 Rabat - MOROCCO

